

THE FINAL EVOLUTION of KAIZER
– KAIZER XX



KAIZER X EXPERT & EXPRESS



EXPERT & EXPRESS

Huvitz

KAIZER X
with Built-in Drilling System

THE FINAL EVOLUTION of KAIZER
– HUVITZ KAIZER XX

With KAIZER XX's absolute high speed, increase
the efficiency and economic feasibility of your business.



KAIZER's Pride

80 Seconds / Fast Control

80 seconds to cut the frame with KAIZER XX.
KAIZER XX has enhanced not only its high speed, but also its stability.

35 Styles of Step-Bevel / Ultimate Flexibility

KAIZER XX supports almost all of frame types with the combination of 7 different step bevel types and 5 edge styles. For example, full or partial, front or rear, bevel top or grooved top, and etc.

2X Built-in Drill Solution / Synergy Effect

Built-in Driller reduces processing time by 50%.
Enjoy the reduced processing time, greater diversity and advanced scalability.

Smarter Performance / Art of Work

Easier to process any kinds of jobs, including complex and difficult ones.
– Smarter algorithm, automatic recognition, intuitive GUI, compatibility, even with crystal clear polishing.



Perfect Fitting

Fitting Size per Frame Material
Efficient and Powerful Retouch

Processing Data

High-Resolution Tracing
Digital Pattern
Image / CAD File Converter
Digital Scan
Hole Editor

Safe Processing

Hydrophobic Processing Mode
Adaptive Chuck/Block Adapter
Powerful Processing Wheel

Lens Blocking

One-Touch Auto Blocking
Auto-Recognition of Lens Diopter / Marking
Auto-Correction of Prism Error


Lens Processing

Various Bevel Shapes
Auto Bevel Position Mode
Editor of Bevel Processing Position
Editor of Variable Asymmetry
Built-In Driller Processing
Auto Grooving Position Mode
Grooving on Scan & Cut Shape
Step Bevel Editor
Scan & Cut Data Editor

5 Insights / Expert

KAIZER has been developed into an innovative function through the expert opinions and researches of high technology for the past 10 years.


We provide the best quality by adopting 5 representative insights that satisfy experts.








Fashionable Sunglasses

As time goes by, frames and materials have been diversifying. And, lens coating with functionalities is easy to slip. With the upgraded KAIZER in high speed with smart functions, it provides breakthrough in these issues.

30%
faster than HPE-8000X



-  Processing Data
 - High-Resolution 3D Tracing
-  Lens Blocking
 - One-Touch Auto-Blocking
-  Lens Processing
 - Variety of Bevel Processing
 - Bevel Position Editor
-  Safe Processing
 - Block Adapter / Adaptive Chuck
-  Perfect Fitting
 - Saving Fitting Size per Frame Material
 - Retouch Function



Processing Data

High-Resolution 3D Tracing of Four Kinds of Frame Material

- 3D tracing with accuracy of 28,000pts is supportive for predefined frame materials like metal, hard plastic, soft plastic and ultem.
- Tracing in ultem mode automatically restores distorted shape, which can be caused while tracing a super-flexible frame.
- Easily possible to change shapes in various ways by using digital pattern function.




High-Resolution 3D Tracing / Restoring Ultem Shape

Lens Processing

Various Bevel Processing for Sturdy Fitting

- Mini bevel, Asymmetric bevel, Semi U bevel, Step bevel
- More accurate fitting by complex processing combining bevel, rimless and groove



Blunt Bevel / Mini Bevel / Asymmetric Bevel / Semi-U Bevel

Manual Position Editor

- Manual position mode easily selects the most suitable position for its frame.
- For general frames, it's convenient to process by auto position mode.



Manual Position Editor

Lens Blocking

One-Touch Auto Blocking

- Auto recognition of optical center or markings by pushing one button.
- Auto compensation of prism error for progressive lense.
- Faster & more efficient job by auto switching of job side and auto transmission of job data after blocking.
- Semi auto recognition mode can proceed with the progressive lenses without markings.



Progressive Lens Recognition

Safe Processing

Block Adapter/Adaptive Chuck to minimize errors & axis distortion

- Newly upgraded block adapter & adaptive chuck, and clamp minimize lens processing errors, and enable safer processing.
- More powerful roughing wheel and smarter sensor with improved algorithm reduce processing time.
- Hydrophobic lens mode's various processing options, clamp pressure, lens rotation speed can be easily set by one button.

Perfect Fitting

Saving Fitting Size per Frame Material

- KAIZER XX supports 4 kinds of material, metal, hard plastic, soft plastic, ultem. Once you register the size standard per material, it's possible to quick processing & fitting.

High Quality of Retouch Function

- Mini bevel is available by retouching bevel height. Adding inclined cut is possible if necessary.
- Depending on the job, you can add grooving, polishing, safety bevel afterwards. Also, you can minimize fitting time by selecting each step such as finishing, grooving or drilling only.

Rimless Glasses

By OMA file conversion, immediate processing without drilling layout editing is available.

32%
faster than HPE-8000X



Processing Data

- Importing OMA File
- File Converter
- Digital Scan
- Digital Pattern
- Hole Editor



Processing Data

Direct Importing OMA File

- Once inputting DCS (OMA) file such as provided from well-known frame brands into memory card, you can easily import and process immediately without design work.

File Converter to increase Design Creativity

- Inputting users' own special designs into memory card, file converter automatically imports frame shapes, drilling information from standard image files, or DXF or CAD files.

Immediate Digital Scan

- Holes and frame shapes can be recognized by high-resolution camera.
- Putting lenses on the lens support, scanning starts without distortion caused by gravity and accurate size of frame shape can be acquired.
- In case of Scan & Cut with drilling, data editor automatically carries out, can delete data for the position that users don't want drill cutting.

Improved Digital Pattern


- With the digital pattern, user can create unique frame shape with the image background of scanned demo lens.
- Users can choose curve or line editing, and apply to any positions.
- Stretching to secure space of near portion in progressive or multi-focal lenses.
- Shaping corner of angulated shape.
- Restoring damaged demo lenses.
- Rotating shape to change position of horizontal axis.
- Stretching shape to specific direction or enlarging and reducing the total size.

Hole Editor to create high-quality glasses


- In hole editing mode for rimless, users can work with the image background of scanned demo lens.
- Easily adding or correcting hole, slot, notch is possible. It supports square and circle shapes as well.
- Users can choose drilling position based on center/edge/box edge, or the distance from a specific hole.
- Users can import and use the saved drilling information. Registering drilling patterns on preset, more efficient editing job is possible.
- By easy click mode, processing for Chemistrie Clip can be operated easily.




Direct DCS(OMA) Import



File Converter / Scan & Cut Shape



Digital Pattern



Hole Editor




Built-in Drill / Chemistrie Clip

Trendy Premium Glasses

KAIZER XX has been upgraded with Improved retouch, grooving functions so it can perfectly create lenses for trendy premium glasses.



40%
faster than HPE-8000X



Lens Processing

- Various Shapes of Complex Processing
- Partial Groove Editor
- Scan & Cut Shape Grooving

Safe Processing

- Powerful Roughing Wheel / Adaptive Chuck

Perfect Fitting

- Retouch Function

Lens Processing

Complex Processing to Perfect Fitting for Various Shapes

- By supporting 4 kinds of complex processing combining bevel, rimless and grooving, users can perfectly finish fitting for unique and luxurious frames.
- Partial groove (grooving + rimless)
- Dual groove (basic grooving + special grooving)
- Hybrid groove (grooving + bevel)
- Partial bevel (rimless + bevel)

Partial Grooving Editor


- Users can designate grooving area & safety beveling. Also, can use and angle guide with zoom in/out for precise edit.

Grooving for Scan & Cut shape


- For unique designs from house brand and handmade frames, KAIZER XX supports grooving of concave shape.



Hybrid Grooving / Dual Grooving / Partial Grooving / Partial Bevel



Setting Groove on Scan & Cut Shape



Adaptive Clamp Chuck

RPA/ RPGA Wheel Types

Safe Processing

Upgraded Roughing Wheel / Adaptive Chuck

- The upgraded roughing wheel & adaptive chuck support powerful and speedy processing with hydrophobic mode.
- RPGA type : Roughing wheel + Hybrid finishing wheel
- RPA type : Wide roughing wheel + Hybrid finishing wheel

Perfect Fitting

Improved Retouch Function

- Improved retouch function enables more efficient fitting for complex processing using grooving or bevel.
- The function dramatically reduces processing time by specifically adjusting the range as well as size and grooving.




Partial Grooving Retouch




Sports Goggles

KAIZER with built-in drill & twice faster speed, customize processing for high curve sports goggles, boutique premium design.



45%
faster than HPE-8000X




-  Processing Data
 - Digital Scan
-  Lens Processing
 - Step Bevel Editor
 - Scan & Cut Data Editor
 - Built-in Drill Processing
-  Perfect Fitting
 - Retouch Function

Processing Data

Digital Scan for High Curve Goggle

- Recognizing step bevel line properly and compensating distortion caused by gravity are the most important jobs in high level goggle processing.
- The digital scan using special lens support can accurately scan frame shapes.
- Step bevel line can be indicated by manual or already processed step bevel lenses.



Demo Lens Digital Scan

Lens Processing

Step Bevel Editor

- 7 types of step bevel are available. (Full step, partial step can be processed on front, rear and both sides.)
- Step bevel edge can be finished in flat, bevel, grooving, blunt bevel and inclined cut.



Step Bevel Editor

Scan & Cut Data Editor

- Using Scan & Cut Editor, users can easily reduce the number of the vent holes in sports goggle design.



Scan & Cut Data Editor

Step Bevel Processing

Built-in Drill

- The high-speed built-in drill enables 2X faster processing for Scan & Cut shape. One-stop processing from drilling to step bevel is possible.


Perfect Fitting

Retouch Mode

- Retouch mode even enhances the efficiency by allowing wheel-edging, drilling or step beveling selectively in smart way.



Retouch Mode



Step Bevel / T-Step Bevel / Beveled Top / Grooved Top

Kaizer Series

Huvitz edging system offers various solutions meeting customers' needs.

HPE-910, Upgraded with Processing Tools & Motor Engine

- Built-in drill with 2X longer durability, reduced processing time more than 50%
- Perfect fitting by step bevel with 7 types and 5 edge styles.
- With more powerful wheel, more flexible adaptive chuck, easy to process high curve and various coating lenses
- Fast and easy for complex processing include variable asymmetric which is adjustable for height.
(Semi-U, Asymmetric, Partial Grooving/Beveling, Hybrid Grooving, Dual Grooving)
- Concave shape grooving, adjusting frame bevel's height, Improved retouch.



HDM-8000, Edging & Drilling at the same time

- Stand alone drilling system which can process rimless frame at the same time with edging.
(HDM-8000 is compatible with 910X, 910N types.)
- Manual angle mode/Scan & Cut mode to process various designs including sports goggles.
- Easy click mode for Chemistrie Clip.



HAB-8000X with Faster Speed & High Efficiency

- By advanced internal lensmeter & camera, it automatically scans sports goggles, demo lenses, tinted lenses without error.
- Auto recognition of optical center & blocking position for multi-focal & progressive lens.
- Auto processing for error compensation for prism distortion and input diopter/axis.
- Accurate recognition of frame shape including step bevel line and hole/slot of rimless frame.
Shape of damaged goggle can be restored as well.
- Smart convenient function including smarter algorithm, adding type for ultem frame, standing tracing and layout in a line.
- Powerful digital pattern & designer function such as creating own design, free editing.




Kaizer Series

Huvitz edging system offers various solutions meeting customers' needs.



Kaizer Configuration

Choose any configuration, Huvitz offers high-quality performance.



HFR-8000X with 3D Scanning includes Faster Speed & Improved Accuracy

- HFR-8000X provides speedy scanning with high-resolution of 16,000pts by adopting the latest motor control platform and digital signal processing technology.
- It supports 3D size compensation function so that can easily recognize concave lenses.
- 3D smart scan without distortion for angle edge, narrow frame and ultem frame.
- By upgraded algorithm, able to trace gripper bar's position, and expect real-time frame shape.

HMB-8000 with Compact Design & High Efficiency LED Light Source

- Adopting high efficiency LED light source which can easily recognize dark-tinted lenses, and adjust brightness.
- Compact design, intuitive graphic, auto power on/off function.

Specification by Each Type

| Product Type | Wheel Type | Asymmetric | Built-in Drill | Step Bevel | Scan & Cut | Chemistrie |
|--------------|-------------|------------|----------------|------------|--------------|------------|
| HPE-910(XD) | 110V RPA* | O | O | O | O (Built-in) | |
| HPE-910(X) | 110V RPA* | O | X | O | O (HDM-8000) | |
| HPE-910(N) | 110V RPA* | O | X | X | O (HDM-8000) | |
| HPE-910(XD) | 220V RPGA** | O | O | O | O (Built-in) | |
| HPE-910(X) | 220V RPGA** | O | X | O | O (HDM-8000) | |
| HPE-910(N) | 220V RPGA** | O | X | X | O (HDM-8000) | |

* RPA type consists of high-speed super wide roughing wheel and hybrid finishing wheel.


** RPGA type consists of high-speed roughing wheel positioned left-hand side of the glass wheel and hybrid finishing wheel.

Network Solutions & Specifications

Huvitz edging system provides optimizing network solutions considering user friendly settings, system configurations and device maintenance.

HERA Intraworks (STD) Optical Store Solution


- Standard version serves as a data server among Huvitz finishing instruments and provides connections up to the capability of PC.
- Its simple interface allows data edition, job assignment, and job import/export even in OMA format which can be automated by options.



HERA Intraworks (RT) Remote Tracing Solution

- RT version is to help establish the simplest remote tracing environment by the acquisition of tracing data and automatic export of it in OMA format for the customer to upload it on LAB-ordering sites.

*Remote Lab receives job data via a Lab Management Software, or can make an economic solution with the combination of email and HERA Intraworks(STD).




HERA Intraworks (HOST) Third Party Compatibility Solution

- Host version based on DCS (OMA) supports data connection between Huvitz devices and third party's devices.

- Host version provides simple third party compatible solution without expensive dcs host solution.

*To be compatible with Huvitz system, 3rd party instruments should support OMA protocol and may require a compatible chuck adaptor.



HPE-910 (Edger)

| | |
|---------------------|--|
| Lens Material | Plastic, Polycarbonate, High index plastic, Glass, Trivex |
| Edging Type | Bevel, Groove, Rimless, Mini bevel, Partial bevel, Partial groove, Dual groove, Hybrid groove, Asymmetric bevel, Semi-U bevel, Step bevel* |
| Edging Position | Front %, Front mm, Rear mm, Base curve, Auto, Manual (with 2D simulation and 3D preview) |
| Edging Options | Polishing, Safety bevel, Safety mode |
| Built-in Drilling** | Hole, Slot, Notch, Scan & Cut shape, Chemistrie Clip |
| Functions | Job manager, Digital pattern, Hole editor, Asymmetric bevel editor, Scan & Cut shape editor, Step bevel editor, OMA import/export, CAD export, DCS/OMA compatibility, Retouch mode (Normal, Drill-only, All) |
| Utilities | LCD tilting Automatic edging room door Edging room illumination SD card storage (Memory card included) Barcode reader interface, Vacuum interface |
| Edging Size | Max. 90mm Min. flat edging : 18.5mm (without safety bevel) Min. bevel edging : 20mm (without safety bevel) |
| Display | 9.7" color TFT LCD(1024x768) with touch screen |
| Dimensions / Weight | 540(W) x 472(D) x 580(H)mm 51.2kg (110V), 50.2kg (220V) |
| Power Supply | AC 100-120V / AC 220-240V 50/60Hz |
| Power Consumption | 1400VA |

* 'N' type doesn't support step beveling. ** 'X' and 'N' types don't support built-in drilling.

| Wheel / Edger Type | RPGA (220V) | | | RPA (110V) | | |
|--------------------------|--|---|---|------------|---|---|
| | XD | X | N | XD | X | N |
| Wheel Combination | High-speed plastic roughing (220V: normal, 110V: super wide) | | | | | |
| | Hybrid finishing (Finishing + Asymmetric) | | | | | |
| | Polishing (Rimless only) | | | Polishing | | |
| | Glass roughing | | | None | | |
| Lens Material | Plastic | | | | | |
| | Polycarbonate | | | | | |
| | Hi-index | | | | | |
| | Trivex | | | | | |
| | Glass | | | None | | |
| Mini Beveling | 0 | 0 | 0 | 0 | 0 | 0 |
| Asymmetric / Semi-U | 0 | 0 | 0 | 0 | 0 | 0 |
| Partial Grooving* | 0 | 0 | 0 | 0 | 0 | 0 |
| Hybrid Grooving* | 0 | 0 | 0 | 0 | 0 | 0 |
| Built-in Drilling | 0 | X | X | 0 | X | X |
| Drilling Machine Support | X | 0 | 0 | X | 0 | 0 |
| Step Beveling** | 0 | 0 | X | 0 | 0 | X |

*No grooving support for glass material. **No step beveling support for glass and plastic, or CR39.

HDM-8000 (Drilling Machine)

| | |
|------------------------|--------------------------------------|
| Hole Type | Hole, Slot, Notch |
| Hole Size | ø1.00~5.00mm |
| Tilting Scope | Automatic, Manual (0~30°) |
| Hole Depth | Max. 6.0mm (0.0mm = through hole) |
| Range of Hole Drilling | ø32.0~75.0mm from lens rotation axis |
| Slot Width | 1.0mm~5.00mm |
| Slot Length | Max. 20.00mm |
| Dimensions / Weight | 193(W) x 483(D) x 342(H)mm / 14kg |
| Power Supply | AC 100-240V 50/60Hz |
| Power Consumption | 100W |

HAB-8000X (Auto Blocker)

| | |
|-----------------------|---|
| Tracing Type | (Refer to the specification of the frame reader) |
| Tracing Mode | |
| Tracing Size | |
| Frame Material | |
| Frame Data Processing | |
| Measurement | SPH: -10D~+10D, Cyl: ±6D (min. step 0.01D) |
| Blocking Tolerance | -0.5~+0.5mm, Axis ±1° |
| Blocking Method | Automatic free location |
| Blocking Pressure | 3.0kgf |
| Lens Recognition Mode | Intelligent, Auto, Manual, Semi-auto (for laser mark) |
| Lens Type Recognition | Single, Bi-focal, Progressive, 3-dot |
| Digital Scan | Shape, Hole/Slot/Notch, Step bevel line |
| Blocking Assistance | Camera zoom, PD/OH live adjustment, Lens shortage check, Block attachment check, Auto brightness control |
| Layout Factors | FPD/DBL, PD (Binocular, Monocular), Cylinder axis, OH (ΔY, Mixed height, Box height), Centering method (Box center, Optical center) |
| Edging Factors | Lens material (Refer to the edger) Edging type (Refer to the edger) Edging position (Refer to the edger) Edging option: Polishing, Safety beveling, Safety mode |
| Functions | Job manager, Digital pattern, Hole editor, Asymmetric bevel editor, Scan & Cut shape editor, Step bevel editor, Image/CAD file converter, OMA import/export, CAD file export, Auto job save & transmission, R/L shape correction, DCS/OMA compatibility |
| Utilities | Tiltable LCD screen SD card storage (Memory card included) Barcode reader interface |
| Display | 10.4" color TFT LCD(1024x768) with touch screen |
| Dimensions / Weight | 300(W) x 470(D) x 560(H)mm / 23kg |
| Power Supply | AC 100-240V 50/60Hz |
| Power Consumption | 75W |

HFR-8000X (Frame Reader)

| | |
|---------------------|--|
| Tracing Type | Automatic 3D binocular tracing |
| Tracing Mode | Auto, Semi-auto (Manual tip-positioning) |
| Tracing Size | Frame ø16.0~92.0mm, Pattern ø16.0~84.0mm |
| Frame Material | Metal, Hard plastic, Soft plastic, Ultem |
| Data Processing | FPD, Frame curve, Circumference, 3D angle, Concave shape (Demo lens) |
| Dimensions / Weight | 248(W) x 320(D) x 190(H)mm / 8kg |
| Power Supply | AC 100-240V 50/60Hz |
| Power Consumption | 32W |

HMB-8000 (Manual Blocker)

| | |
|---------------------|---|
| illumination | White LED source Light intensity adjustment Automatic power-saving mode |
| Dimensions / Weight | 177(W) x 184(D) x 206(H)mm / 2kg |
| Power Supply | 5V DC 3.5A |
| Power Consumption | 2.5W |

Designs and details can be changed without prior notice for the purposes of improvement.