



DAUD JEE MFG.CO

We are Part of Saving Vision



Ophthalmic
Instruments



www.daudjeemfg.com

Catalog Number Prefix	Product Group	Page Number
D10	Speculums & Retractors	04-11
D20	Knives, Spuds, Blade Breakers,Punches & Handles	12-14
D30	Lasik/Lasek/Dalk,Curettes,Spatulas,Choppers, Manipulators,Loops,Spoons,Hooks,Fixation Rings, Muscle Hooks,Corneal Markers & Calipers.	16-48
D40	Scissors	50- 60
D50	Forceps	62- 100
D60	Needle Holders	102- 106
D70	Lacrimal Instruments, Cannulas, Quickswitch I/A System Infusion Handles Bimanual Hand Pieces, Rongeurs, Nasal Speculums, Elevators, Chisels, Mallets & Probes.	107- 113
D80	Bipolar Forceps.	114
D90	Sterilization Cases & Trays, Diagnostic Instruments	115-118


D10
● Speculums
● Retractors

D20
● Knives
● Spuds
● Blade Breakers
● Punches
● Handles

Speculums

Barraquer Wire Speculums


D10-5000 standard wire 0.8mm, 10mm blades
D10-5009 heavy wire 1mm, 15mm blades
D10-5010 standard wire 0.8mm, 15mm blades



Barraquer Wire Speculums


Wire Speculum
D10-5355 10mm blades (infant)
D10-5356 15mm blades

Grandon-Barraquer
D10-5011 15mm blades



Temporal Approach


Wire Speculum
with spring on nasal side
D10-5012 15mm closed blades
D10-5015 15mm open blades



Kratz-Barraquer


Wire Speculum
with open blades
For unobstructed access to the
cataract incision with phaco or I/A tip.

D10-5049 heavy wire 1mm, 15 mm blades
D10-5050 standard wire 0.8mm, 15mm blades
D10-5051 standard wire 0.8mm, 13mm blades
D10-18005/A standard wire 0.8mm, 15mm blades
D10-18005/B heavy wire 1mm, 15 mm blades



Barraquer Solid Blade

Wire Speculum
D10-5020 Barraquer, 10mm blades
D10-5030 Barraquer, 15mm blades
D10-5031 Feaster, 20mm blades



Speculums

Sauer

Lid Speculum
D10-5300 for infants, 11mm blades
D10-5302 for premature infants, 2mm blades




Frankel

Infant Lid Speculum adjustable
for precise control of
infant eyelid retraction.
D10-5310 10mm blades



Infant Wire Speculum

D10-5320 Blades size 5mm, Open
D10-5325 Blades size 5mm, Close
D10-5330 Blades size 6mm, Solid



Alphonso

Infant Lid Speculum finger grips and
interlacing blades to facilitate insertion
one blade 5mm, one blade 7mm
D10-5340



Infant Lid Speculum


Lid Speculum closed wire blades
gently curved spring
D10-5350 4mm blades



Speculums

Murdoch

Eye Speculum with locking screw
D10-5406 #2, 12mm blades
D10-5407 #3, 15mm blades



Mellinger

Eye Speculum self retaining
for unobstructed access to the
cataract incision with phaco or I/A tip.
D10-5411 15mm blades



Williams

Eye Speculum with locking screw
D10-5500 11mm blades
D10-5510 15mm blades



Knapp

Eye Speculum
D10-5511 13mm blades



Lancaster Speculum

Eye Speculum
adjustable mechanism
with locking nut
D10-5600 16mm blades



Speculums

Lester - Burch Speculum

Eye Speculum with flanges on
blades for retaining eye lashes
adjustable mechanism
with locking screw
D10-5620



Speculums


Lieberman Speculum

Eye Speculum
Adjustable mechanism
non-reflective finish
D10-5666 5mm blades, Infant
D10-5677 10mm blades, Small
D10-5678 15mm blades, Standard



Jaffe

Wire Lid Retractor
D10-8000 pair, 15mm blades (large)
D10-8001 pair, 13mm blades (small)




Castroviejo Speculum

Eye Speculum
adjustable mechanism
D10-5700 5mm blades
D10-5701 10mm blades
D10-5702 14mm blades
D10-5703 16mm blades



Bangerter Speculum

D10-1230 15mm blades
D10-1240 10mm blades



Stevenson

Lacrimal Sac Retractor
3x3 blunt prongs, adjustable
D10-9500



Stevenson

Lacrimal Sac Retractor
solid blades with serrated
edge adjustable
D10-9520 13mm blades



Retractors


Schocket

Scleral Depressor double-end
with pocket clip 6mm wide bar
and 2.5mm tear drop
D10-7825



Josephberg-Besser

D10-7822



Schepens

Orbital Retractor
D10-7880



Senn

Senn Miller Lid Retractor
D10-8100
Sharp, 7x8mm, 16 cm 6¼"
D10-8120
Blunt, 7x8mm, 16 cm 6¼"




Senn

Senn Miller Lid Retractor,
Double Ended
D10-8150 Sharp
D10-8160 Blunt



Desmarres

Lid Retractor
thin solid blades
D10-8200
size 0, 11mm wide
D10-8210
size 1, 13mm wide
D10-8220
size 2, 15mm wide
D10-8230
size 3, 17mm wide



Retractors

Jaeger

Lid Plate
D10-8500 stainless steel, Small
D10-8510 stainless steel, Large



Wright

</

Knives & Spuds

Tooke

Corneal Knife
straight blade
with curved cutting edge
D20-3650



Gill

Corneal Knife
curved cutting edge
D20-3700 straight
D20-3710 light curve



Fukasaku

"Hockey" Knife
sharp 3mm front edge
semi-sharp 6mm
bottom edge
D20-3735



Foreign body spuds

D20-4000
Golf Club, angled,
lightly rounded tip
D20-4050
Francis, angled,
flag shaped tip
D20-4100 Ellis,
lightly curved and
rounded tip



D20-4000



D20-4050




D20-4100



Spud & Gouge

in reversible handle for
removing foreign bodies.
D20-4600



Handles/Blade Breakers

Castroviejo

Blade breaker and holder
with sliding lock.

D20-5000



Swiss Model

Blade breaker and holder
spherical tip.

D20-5150



Chuck Handle

round, serrated

D20-1910 10cm

D20-1920 7.5cm



Bard-Parker Handle

for disposable slot blades

D20-1950 #3 with gauge

D20-1960 #3



Barron Handle # 3

D20-1971



D30

Curettes

Spatulas/Choppers/

Manipulators/Loops

Loops/Spoons/Manipulators

Hooks/Manipulators

Fixation Rings

Muscle Hooks

Corneal Markers


Calipers

Curettes/Spatulas

Meyerhoefer

Chalazion Curette

D30-1000 size 0, 1.5mm dia.
D30-1010 size 1, 1.75mm dia.
D30-1020 size 2, 2.25mm dia.
D30-1030 size 3, 2.5mm dia.
D30-1040 size 4, 3.5mm dia.



Double End Currrette

D30-1041 size 1.5mm / 2.5mm dia.
D30-1042 size 2mm / 3mm dia.
D30-1043 size 2.5mm / 3.5mm dia.



Capsule

Polishing Curette

1.75mm diameter open ring
D30-1100 Nightingale, semi-sharp, posterior edge
D30-1102 Shepherd, semi-sharp, anterior and posterior edges



Culler

Iris Spatula

D30-2100 1mm wide
D30-2110 2mm wide
D30-2111 nettleship iris repositor



Wheeler

Double-End Spatula
lightly curved

D30-2120




Barraquer

I

Endothelial Stripper

1mm wide, thin beveled tip
D30-1830
smooth surface
D30-1835
sosenwasser with textured surface



Sloane

LASEK Flap Repositor
gently curved spatulated tip
with beveled smooth edges
for repositing and smoothing
the epithelial flap after laser
ablation
D30-1845



Camellin style

LASEK Trephine &
alcohol retaining well
D20-7830
8mm trephine and
8.5mm well




Camellin Style

LASEK Alcohol Well
For retaining alcohol solution
on the corneal epithelium.
D30-1820 8.5mm diameter
DT30-1820 8.5mm diameter titanium
D30-1822 9.5mm diameter
DT30-1822 9.5mm diameter titanium



Anis

Capsule Scrubber flexible
wire textured disc
D30-1111




Moran

LASIK Enhancement
Spatula Fukasaku spatula
and a modified paton spatula
D30-252

Lasik/Lasek/Dalk

Hersh

Flap spatula beveled tips
a combination of a 6mm
spatula and 10mm spatula
D30-2526



Binder

Flap elevator 1mm wide 12mm
semi-sharp edges
D30-2528



Fukasaku

LASIK Spatula rounded anterior
and flat posterior surface with
beveled tip for lifting and
dissecting the corneal flap.
D30-2530




Guimaraes

LASIK spatula
gently curved, bi-convex
semi-sharp edges
D30-2532



Katzen

LASIK Flap "Unzipper"
combines a modified
sinskey hook and a
bi-convex spatula
D30-2534



Maddox

LASIK Spatula
combines a flat spatula and
a cylindrical spatula for
elevating and repositing the flap
D30-2535



Spatulas /Choppers


Phaco chopper

Chopper 10mm angle to tip
blunt cutting tip
DT30-2180.8
Chopper 10mm angle to tip
rounded front, 90° cutting edge.
DT30-2180.25
Delicate 10mm angle to tip
short tip, 90° cutting edge.
DT30-2180.29




Mics Phaco Pre-chopper

For bimanual nucleus pre-chopping
D30-2323 Alio-Rosen,
with two defined edges.
D30-2324 Alio-Scimitar, curved
with blunt tip and sharp inferior edge.



Milder

Phaco Chopper
60° offset inferior wedge
for horizontal chopping and a
beveled distal tip for vertical chopping
D30-2343 for use in left hand
DT30-2343 for use in left hand titanium



Brierley

Nucleus Splitter
with incision sealing shaft
D30-2380
DT30-2380 titanium




Dodick - Kammann

Nucleus Chopper
for bimanual nucleus division.
D30-2382 1.5mm long tip, str
D30-2383 1.5mm long tip, angled
<

Spatulas/Nucleus Manipulators


Gimbel

Nucleus Spatula
round, smooth
angled 30° at
10mm from the tip
0.25mm diameter
D30-2312




Knolle

Nucleus Spatula
malleable stainless steel,
flat 0.5mm wide
D30-2350



Koch

Nucleus Spatula
gently curved duckbill
shaped tip
D30-2354 smooth
D30-2355 with notches





Escaf

Nucleus Manipulator
flattened olive-shaped tip angled at 45°
D30-2347
DT30-2347 titanium

Akahoshi

0.5mm diameter bulbous tip,
tip length 1.5mm 45° angled
shaft, tip to angle length 10.0mm
D30-2357
DT30-2357 titanium



Escaf

Nucleus Sustainer
smooth ball tip textured inferior surface
D30-2359

Lieberman

MicroFinger Nucleus Manipulator
D30-2360 for use in right hand
D30-2361 for use in left hand



Spatulas/ Manipulators

Lippman

Nucleus Spatula
0.5mm diameter round
shaft with blunt forked tip
D30-2430



Girard

Synechiaie Spatula
0.5mm wide with
semi-sharp tip
D30-2500



Obstbaum

Synechiaie Spatula
0.25mm wide with
semi-sharp tip
D30-2520



Benedetto

Spatula 10mm long,
smoothly polished
flattened spatula
with beveled notches
D30-2522



Elsching

Cyclodialysis Spatula
flat, 1mm wide.
D30-3000



Spatulas/ Manipulators

Castroviejo

Synechiaie Spatula
0.5mm wide with
semi-sharp tips.
D30-3150



Lens Loops

Keener

Lens Divider
with three preformed snare
wire loops
D30-3900 complete
D30-3910 wire loops



Lens Loop

D30-4000
Lewis, small
D30-4010
Lewis, Large
D30-4030
New Orleans
D30-4100
Wilder, serrated



Kirby

Lens Loop and
Muscle Hook
D30-4105



Morrison

Lens Loop and
Probe
D30-4120



Lens Loop

serrated, 3mm wide
D30-4140
Kansas, light curve
D30-4141
Keener-Arlt,
strong curve



Lens Loops/Spoons

Alfonso

Cutting Platform
Used as a cutting
block for

Spoons/Hooks

Wells

Enucleation
Spoon
D30-4450




Graether

Collar Button
micro iris retractor
D30-4900 straight,
with guard
D30-4910 angled



Keuch

Small Pupil Dilator
with guard for enlarging
a small pupil through the
sideport and primary incision.
D30-4950




ICL Manipulator

D30-5011



Bonn

Micro Iris Hook
extra delicate blunt,
round handle with guard.
D30-5050



Hooks/Fixation Rings


Fixation Hook

double, sharp
D30-6000
D30-6010




Twist Hook

for scleral fixation
D30-6100 right
D30-6110 left




Silverman

Fixation Pick
blunt cone shaped tips
spaced 3mm apart
for atraumatic scleral fixation
D30-6130



Globe Ring

Globe Fixation Ring
with concentric grooves
D30-6153 for refractive procedures
D30-6154 for cataract procedures



Lu - Mendez

LRI Guide and Fixation Ring
with delicate blunt teeth marks 10°
increments 11mm I.D.
D30-6158
DT30-6158 titanium



Fixation Rings

Gimbel

Stabilization Ring
with swivel handle
13mm diameter with blunt teeth
on both sides
For stabilizing either the right or
left eye during LASIK.

Muscle Hooks


Stevens
Tenotomy Hook
blunt, small
D30-6600



Helveston
"Teaser" Hook
extra delicate, blunt
D30-6610



Graefe
Muscle Hook
D30-6710 size 1
D30-6720 size 2
D30-6730 size 3



Gass
Retinal Detachment
Hook with oval hole
D30-6740



Helveston
Muscle Hook
extra delicate,
with blunt spatulated tip
D30-6760 large
D30-6761 small



Muscle Hooks

Helveston
Muscle Hook
extra delicate, with blunt
spatulated tip
D30-6778 8mm
D30-6780 10mm
D30-6782 12mm



Jameson
Muscle Hook
D30-6800 small
D30-6810 large



Corneal Markers

Mendez

Degree Gauge 12mm I.D.
marked in 10° increments
for accurate axial
alignment on the cornea
DT30-7900 titanium
D30-7900




Bores

Axis Marker
DT30-7910 titanium
D30-7910




Henderson

Degree Gauge 11.7mm I.D.
Beveled face with markings
in 5° increments 4 notches
on internal edge to align
with cardinal positions
DT30-7904 titanium




Henderson

Toric IOL Marker
two 2.5mm radial
marks extending to 11mm
(used with ST30-7904)
DT30-7912 titanium



Henderson

Henderson Alignment Marker
with atraumatic disc-shaped
marking surfaces
at 3, 6 & 9 o'clock
DT30-7908 titanium



Corneal Markers




PRECISE, CONTROLLED DELIVERY, TIME AFTER TIME

Y-TIP PLUNGER ROD


- Rounded, Capsular-friendly dimensions assist lens manipulation post-implantation
- Distinctive blue color designed for easy visibility in the cartridge and reduced microscope reflections

DT50-7791



CARTRIDGE LOCK

Quick and easy snap-in design, for secure cartridge loading, assures correct position of cartridge every time.



SCREW-STYLE

Consistent insertion forceps result in predictable lens delivery from start to finish.

High quality titanium material and durable construction for safe, reliable implantation

D40

- Scissors

Metzenbaum

Dissecting Scissors
gently curved blades
blunt tips
5 3/4" (14.5cm)
long

D40-9040

DTC40-9040 Tungsten Carbide

DSC40-9040 Super Cut



Utility Scissors

one pointed and one
blunt tip
5 3/4" (14.5cm) long
straight

D40-9060

DTC40-9060 Tungsten Carbide

DSC40-9060 Super Cut



Enucleation Scissors

blunt tips
D40-9910
medium curve
D40-9920
strong curve



Black Coated Super Cut Scissors

black coated super cut Iris scissors
DSCB40-7400 Str
DSCB40-7410 Cvd



D50

•Forceps

D60

- Needle Holders

D70

- Lacrimal Instruments
- Bimanual I/A
- Quickswitch I/A System
- Infusion Handles
- Roungeurs
- Nasal Speculums
- Elevators
- Chisels
- Mallets
- Probes

