

SLIT LAMP (3 STEPS)

O-SL-11

Diopter Adjustment Range :	-6D to +6D
Total Magnification :	10x, 16x, 25x
Filters :	Cobalt Blue, Green, Yellow, Red
Slit Rotation :	0-180
Eye Piece Magnification :	12.5x
Illumination Type :	Halogen
Slit Width & Slit Length :	1-12mm
Slit Aperture :	0.2mm, 01mm, 03mm, 06mm, 09mm, 12mm
IPD Range :	45 - 75 mm
Illumination :	12V 50 W Halogen
Power Supply :	220V (50-60Hz)
Weight :	20Kgs

HIGH QUALITY
MEDICAL EQUIPMENTS BY

