

Product 1

**Reusable
Video**

Laryngoscope

Item: HYHJ-KC

HYHJ series

Video laryngoscope

Application:

Anesthesia
ICU
Emergency Center
Pediatrics
Ambulance
Intubation training
.....



Video laryngoscope

Power Button

ON/OFF&freeze

Display

3 inch HD
screen,UP/DOWN,LEFT/RIG
HT turning

Alloy handle

Ergonomic design,easy handle

Medical stainless steel

suitable for different
kind of disinfection.reusable



USB interface

charging and data transfer



Shooting

short press for shooting,
long press for recording


2 million pixels
CMOS

HYHJ-KC is Reusable video laryngoscope with 5 blades and 1 display. It can separate the display and replace different size blades.

Reusable Video Laryngoscope: The Tip Blade Size of HYHJ-KC

Tip Blade Size	length (mm)	height (mm)	Width (mm)	Max width of insert part (mm)	Mini width of insert part (mm)	Ref. Picture	Remark
Tip Blade 0 (L Adult)	153	54.5	29.23	22.41	14.37		MAC 4
Tip Blade 1 (M Adult)	140	49	29.23	21.37	14.46		MAC 3
Tip Blade 2 (S Adult or L Child)	122	34	29.23	21.4	13.56		MAC 2
Tip Blade 3 (M Child)	77.4	22.74	18.35	14.5	10.34		Miller 1
Tip Blade 4 (S Child & Infat)	59.9	17.46	18.55	14.7	10.6		Miller 0

Remark: Tolerance: $\pm 5\%$.





HYHJ-KC
complete spare parts

Suitable for you
Adapter and Plug

Shared Charger 110-240 v



Chinese Plug



Europe Plug



UK Plug





HYHJ-KC Packing Case

