

# AVENUE SYNC

---

Easy solution to support Aurora IQ handheld fundus camera integration to PACS and worklist server.

OPTOMED CORPORATION

# Avenue Sync features




Avenue Sync integrates Aurora IQ camera seamlessly with hospital systems and enables more efficient work with existing tools.



Avenue Sync transfers:

- DICOM worklist from hospital information system to Aurora IQ camera
- DICOM images from Aurora IQ camera to hospital image archive.



Worklist

Image



DICOM Worklist

DICOM Image

# Extra features



Tools for adding, editing and deleting study information



Supports different workflows



Possibility to view images



Simple UI for effective use