

OPTOMED SOFTWARE

AVENUE

AVENUE AI SERVICE TRAINING DECK

Copy approved for external use

Avenue AI Service


Avenue AI is our gateway to artificial intelligence services. With the cloud-based system, we are able to integrate different AI algorithm providers' services with our cameras and screening software.

- create and review agreements
- create user account
- register Aurora IQ to the account
- follow statistics of Avenue AI usage.

AVENUE

Registration webpage

- Go to Optomed Avenue Service web page
- Select Register



AVENUE

Terms and conditions


- Read Terms and conditions
- Select Agree to terms
- Select Continue



AVENUE

User details

- Fill text fields
- Password must be at least 8 characters
- Password must contain upper, lower and number characters
- Select Next



Optomed
AVENUE

Create an account

Fill out the user profile information. Fields with * mark are required information.


Email	<input type="text" value="john.smith@optomed.com"/>	Address	<input type="text" value="100 Main St"/>
Name	<input type="text" value="John Smith"/>	City	<input type="text" value="Dubuque"/>
Password	<input type="password" value=""/>	Country	<input type="text" value="Ireland"/>
Re Password	<input type="password" value=""/>	Phone	<input type="text" value="+353 20 742 3333"/>

AVENUE

Organization details


- Add Organization name and VAT Number
- Select Register

AVENUE



Confirmation


- Account created successfully
- Wait until Optomed personnel accepts your registration



AVENUE

Login

- Go to front page
- Add username and password
- Select Login



AVENUE

Profile settings


Profile settings

- Change password
- Update and save account information
- View Agreements
- Register Cameras and products

Left side column


- Edit users for administrators
- Help center
- Log out

AVENUE



Register camera

- Select Add camera or product
- Fill fields
 - Product code – Aurora IQ serial number
 - Product name
 - Product description
- Select Register new camera




AVENUE

Register camera

Select camera


- User can add multiple products but only one can be in use at the time



AVENUE

Agreements


- Accepted agreements can be viewed and downloaded
- New terms and agreements can be accepted



AVENUE

Help center

- Contact details
- Avenue service details



AVENUE

OPT**MED**
CORPORATION