

assistTrend

The perfect assistant's stool.



OPHTHALMIC
SURGERY



AESTHETIC
SURGERY



DENTAL
SURGERY



ENT
SURGERY



MICRO
SURGERY



ORAL &
MAXILLOFACIAL

Order No. 25000000


ADVANTAGES



- High-grade, comfortable and easy-care upholstery fabric to meet even the most stringent hygiene requirements
- Sterile, comfortable height adjustment via foot-triggered gas pressure spring
- Excellent movement characteristics thanks to double casters and 5-arm, chrome-plated aluminum chassis
- Reduces typical work-related health problems: muscular tension in the neck and lumbar areas
- Our unique dynamic seating system supports correct sitting postures while protecting your back



Comfortable height adjustment
via foot triggering.



Optional: assistTrend with armrests.

No. 33000000



SPECIFICATIONS

Chassis		Max. Load	110 kg
- Five legged chassis	ø 600 mm		
- Wheels	ø 50 mm	Seat/backrest	
Lifting Column		- Seat/Back adjusting	-6°/+3°
- Height adjustment	54–73 cm	Upholstery	Stamskin Top
- For body height of	1.75 m–1.95 m	Standard color	black
Weight	12 kg	Subject to technical modifications.	

