

## EQUIPMENT



Flexible, adaptable and sturdy – the backrest moves together with the footrest (55°) as the chair is reclined to a Trendelenburg position.



Ergonomic vertical adjustment of the chair back eliminates unnecessary sliding and rubbing. The patient can easily and comfortably lean back while experiencing complete lumbar support.



Stylish colors available for the motionLine. Standard color anthracite 197.



Ergonomically formed, durable, non-sliding rubber grips.

Comfortable, height adjustable neck support pillow as standard version.



Storage basket for patient medical records and belongings.



## SPECIFICATIONS

### Chassis

- Laser frame profile 5 mm
- Chassis length 890 mm
- Chassis width 528 mm
- Wheel diameter  $\varnothing$  300/125 mm

### Lifting Column

- 3 parts support frame
- Total length 1,710 mm
- Width 520 mm
- Upholstery Height 40/60 mm
- Back/foot adjustment 55°

- Armrests 0°/+180°
- Seatrest adjustment to Trendelenburg position 25°

Max. Load 135 kg

Weight 56 kg

Upholstery colors Tundra Skai  
Standard color anthracite  
F6461197

Subject to technical modifications.

