



**Create a Doctor-Patient Dialogue with iNview.
Leverage the power and convenience of the Apple
iPhone with the trusted quality of Volk optics.**

Quickly & effortlessly capture fundus images for visualization & patient education. Helps facilitate patient discussions related to disease progression and treatment plan.

- > Free mobile application available in the Apple App Store (search Volk iNview)
- > Offers 1 Megapixel resolution with a static 50° field of view
- > View the peripheral retina dynamically out to 80°
- > Available auto-capture & forced-capture imaging modes
- > Mydriatic; requires minimum 5mm pupil
- > HIPAA-compliant storage and export from iPhone to PC or Mac
- > Compatible with Apple iPhone SE/6s/6/5s and iPod Touch (6)

