


8 Position Landolt C Near Vision Card (16 inches/40 centimeters)



8 Position Landolt C Near Vision Card (16 inches/40 centimeters)

Rating: Not Rated Yet
[Ask a question about this product](#)

Description

[Click here to view the pdf catalogue](#)

8 Position Landolt C near vision card with C openings at up, down, left, right and each 45° angle. Near vision card with proportionally spaced (logMAR) lines; line sizes range from 20/400 to 20/10 (6/120 to 6/3) equivalent. Attached 16-inch (40 cm) measuring cord ensures proper screening distance. 8 Position Landolt C printed on both sides in scrambled order for repetitive testing. 7" x 9" (17.8 cm x 22.9 cm).
..... 779000