

CHAIN



P/N: 1056

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

- Vision Therapy System for use in a doctor's office and/or for use by the patient at home to help prevent patient regression between office visits
- Three-dimensional Chain image used to strengthen the binocularity system, to provide Base-In & Base-Out training and to provide the patient with an introduction to peripheral binocular experiences
- Empty Chain image allows for objects or another Vectograph to be viewed behind the Chain for additional therapy exercises with increased depth
- Total accommodation range of 40 diopters
- Includes: Protective Therapy Binder with pen, Doctor and Patient Manuals, Therapy Record Form & Standard Polarized Viewer

?

Doctor and Patient Manual

[Printable Doctor PDF](#) [Printable Patient PDF](#)