

MeccanotticaMazza Ophthalmic Unit 10.460.00



MeccanotticaMazza Ophthalmic Unit 10.460.00

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Refraction unit with automatic sliding table for 2 instruments - electric up/down movement of the instruments table and rotation system for patient chair - console with digital control panel with 1 drawer and trial lens tray.

Possibility of adding a second table for 3rd instrument with manual or automatic movement.

Column with adjustable LED lamp and projector support. Patient chair with electric up/down movement.

TECHNICAL DATA

- Power supply 230Vac 50Hz (230Vac 60Hz - 110Vac 60Hz)
- Power consumption 700W
- Auxiliary outputs 230Vac 50Hz 200W (230Vac 60Hz - 110Vac 60Hz)
- Working table voltage 3-12Vac 230Vac 50Hz 200W (230Vac 60Hz - 110Vac 60Hz)