

SLIT LAMP (5 STEPS)

O-SL-21

Diopter Adjustment Range:	-6D to +6D
Total Magnification	: 8x, 12x, 20x, 30x, 40x
Filters	: Cobalt Blue, Green, Yellow, Red
Slit Rotation	: 0-180
Eye Piece Magnification	: 12.5x
Illumination Type	: Halogen
Slit Width & Slit Length	: 1-12mm
Slit Aperture	: 0.2mm, 01mm, 03mm, 06mm, 09mm, 12mm
IPD Range	: 45 - 75 mm
Illumination	: 12V 50 W Halogen
Power Supply	: 220V (50-60Hz)
Weight	: 20 Kgs

HIGH QUALITY
MEDICAL EQUIPMENTS BY

