

Video Laryngoscope with disposable blades

Item: **HYHJ-1320**



Applicable range:

Suitable for anesthesia, emergency department and ICU, pediatrics, ambulance etc, for clinical endotracheal intubation, simulated practice, clinical teaching, difficult airway intubation

Standard 3-inch display

high definition,
strong practicability

Power switch

One button start
with fixed function

Border

medical 304 stainless steel,
high strength, long life

Screen rotation function

Front and rear rotation angle 0-155°,
the left and right rotation angle 0-275°

One-time use of laryngoscopy

different specifications to
meet different requirements,
aseptic packaging, one-time use

Photo recording

one button take pictures
with video function


High-definition camera system

double anti-fog, intuba



HYHJ-1320 is include 3 Tip blade size for each set.





AV output

Data transmission, charging

Reset hole


Digital memory card





HAIYE





HYHJ-1320 Packing



HYHJ-1320 and HYHJ-1330I/II/III are use same size disposible blade (T1/T2/T3) for adults and children.