

Chart Projector CCP-7000



CHAROPS

The CHAROPS CCP-7000 charts with high quality, allow precise acuity test for you.

The CCP-7000 delivers its promises by utilizing the total refraction system to perform any test you want accurately and efficiently. The various kind of chart of the CCP-7000 can be changed more precisely within very short time.



With world wide adaptive charts of the CCP-7000, You can choose any test you want.

Clear and Fast

High Quality Charts

The CCP-7000 has high resolution power of line, 50lines/mm, so provides neat and clean graphics and a fast chart changing speed marking it easy for operators to use and comfortable for customers.



Easy and Convenient


The best networking

Any refraction system including the CHAROPS CDR-7000 can incorporate CCP-7000 to provide fast and efficient vision tests.

Various and Perfect

Wide range of features

The CCP-7000 with 41 charts including Red/Green and Polarizing filters are the new solution for Binocular balance test, Stereo test, Aniseikonia test, Fusion test.



SPECIFICATION

Chart	41 chart 34 mask, Red/Green & Polarization Filters
Projection Distance	2.5~8m
Chart Rotation Speed	Average 0.15 sec
Projection Magnification	30x at 5m
Power Saving	Automatic Switch off (10 min)
Program	2 programs with a maximum of 30 charts each

Tilt Angle	15°
Power Supply	110-120V / 200-240V, 50 / 60Hz
Power Consumption	Max 0.4A
Lamp	6V / 30W (Halogen)
Dimension	367(W) X 219 (D) X273 (H)mm
Weight	6.6kg

Designs and details can be changed without prior notice for improvements.

HUVITZ
SHANGHAI

SHANGHAI HUVITZ Co., Ltd
No.188, Yuanguang RD, Baoshan,
Shanghai, China, 200444
Tel / +86-21-6102-8279
Fax / +86-21-3630-7064
E-mail / alex@huvitz.com

Distributed by

