

CHART PROJECTOR CCP 9000



To perform precise acuity test with high quality charts, chart projector can be the right answer

CCP-9000 delivers its promises by utilizing the total refraction system to perform any test you want accurately and efficiently. The various kind of chart can be changed more precisely within very short time.



Clear and Fast

High Quality Charts

CCP-9000 has high resolution power of line, so provides neat and clean graphics and a fast chart changing speed marking it easy for operators to use and comfortable for customers.

Bright and Efficient

LED light

Semi-permanent lifetime LED light source can provide brighter and clearer charts than conventional bulb light.

Easy and Convenient


Best Networking

Any refraction system including OZMA digital refractor can incorporate CCP-9000 to provide fast and efficient vision tests.

Various and Perfect

Wide range of features

CCP-9000 with 41 charts including Red/Green and polarizing filters are the new solution for Binocular balance test, Stereo test, Aniseikonia test and Fusion test.



SPECIFICATION

Chart	41 chart 34 mask, Red/Green & Polarization Filters
Projection Distance	2.5~8m
Chart Rotation Speed	Average 0.15 sec
Projection Magnification	30x at 5m
Power Saving	Automatic Switch off (10 min)
Program	2 programs with a maximum of 30 charts each

Tilt Angle	15°
Power Supply	110-120V / 200-240V, 50 / 60Hz
Power Consumption	Max 0.6A
Lamp	LED 4W
Dimension	270(W) x 182(D) x 230(H)mm
Weight	3.5kg

Designs and details can be changed without prior notice for improvements.