



Indirect Ophthalmoscope WITH LED WIRELESS

Model : MEDVIC-IOL-1000

Description

- *Compact and Light Weight
- *Has all Pupil Features
- *Stereo Optical System
- *Brilliant illumination
- *Availability of Cobalt Blue and Green Filters
- *Simple & Smooth controls for adjusting the Headband
- *Integrated Teaching Mirrors in Left and Right,

TECHNICAL SPECIFICATION

Type	Rechargeable Led Indirect Ophthalmoscope Wireless
Papillary Distance	45-65 mm
Focus Distance	450mm
Aperture Size	4.0mm, 3.0mm and 1.2mm
Image Size	80.0mm, 60.0mm and 25.0mm
Illumination Field	Good, Clear Circular Field
Lamp	3W Warm White Led
Intensity Control	Continuously Variable
Filters	Yellow, Blue, Green
Power Supply Unit	
Power Input	110V/220V AC
Output	3V DC

Weight	Net 1.5Kg and Gross 3.5Kg
STANDARD ACCESSORIES	
Carrying case	
Charger	
Scleral depressor	
Fundus chart	
Four pencil	
Instruction manual	
OPTIONAL ACCESSORIES AT EXTRA COST	
20 Diopters Aspheric viewing lens	

