

CS-518

OPHTHALMIC UNIT & CHAIR

Space save type: the minimum room width request is 155cm

SPECIFICATIONS

Table Size.....	845mmX405mm
Total Size.....	1000mmX752mmX2100mm
Voltage.....	Main:220V±10% 50Hz; 110V±10% 60Hz
Power.....	Table top: 12V 200W
Net Weight.....	220kgs

FEATURES

Arm

Can be moved up and down or spin around.

Electric lift chair

Ultra-luxury European-style.

Sliding table

Table could slide horizontally, the total length of the machine is only after the table is slid out.

Optometry energy-saving lamps

Elegant, energy-saving.



CE FDA

