



**D903
D903A**

AUTO LENSMETER

COLOR TOUCH SCREEN

FEATURES

- Colorful touch screen, view angle adjustable
- With printer and PD, inner UV part
- Automatic recognition of single and progressive lens
- Strong aluminum case for single packing

SPECIFICATIONS

Measure Scope

- Sphere.....0 to ±25D(0.01D, 0.12D, 0.25D steps)
- Cylinder.....0 to ±9.99D(0.01D, 0.12D, 0.25D steps)
- Axis.....0-180°(1° steps)
- ADD.....0 to 9.99D (0.01D, 0.12D, 0.25D steps)
- Prism.....0 to 15Δ(0.01Δ, 0.12Δ, 0.25Δ steps)

Mode

- Cylinder.....+,+/-,-
- Prism.....X-Y,P-B
- Contact lens.....Soft/hard
- Measuring mode.....single/progressive/
automatic recognition contact lens

Others

- Lens Diameter.....Φ20 to Φ108mm
- PD.....45 to 90mm(0.25mm steps)
- Speed.....0.1S
- ABBE.....30 to 60
- Display.....D903:TFT LCD (5.6")
D903A:TFT LCD (7")
- Printer.....Thermal print
- Dimension.....222(L)x161(W)x372(H)mm
- Package.....430(L)x235(W)x295(H)mm
- Gross weight.....7Kgs
- Power.....AC220V/50HZ; AC110V/60Hz



D903(5.6" touch panel)

CE FDA



D903A(7" touch panel)



Standard Packing: Aluminum Box