


ZN-270 LED VISUAL CHART 27" HD

POLARIZED SCREEN

- Big size monitor, clear display
27" TFT-LED monitor, 1920x1080 pixels, 200CD/m²
- 3D polarization vision chart
Apply to binocular balance test and stereo visual function test
- Optotypes can be changed randomly



CE FDA

Both table stand and wall mount are as standard accessories

SPECIFICATIONS

LED monitor.....27-inch TFT Panel, 1920x1080 pixels, 200CD/m²
 Optotype.....E, C(Landolt),Alphabet, Numbers, Children
 Visual unit.....Decimal,Snellen(ft,20/20),Snellen(m,6/6)
 Dis.....1M~6M (in steps of 0.1m)
 Mask type.....Single mask, Horizontal mask, Vertical mask,
 Red/Green mask, White/Black mask,
 Contrast adjustment, Red/Green adjustment

Special test masks.....Ishihara Test, Cross cylinder,
 Astigmatism Test, Cross Grid,
 Fixation point Test, Red/Green Panel,
 Worth Test, Red/Green Balance,
 Traffic light, Coincidence Test, Stereo Test,
 Schober Test, ETDRS
 Other functions.....Ramdomization Optotypes Playing function,
 R/G setting function, Video playing function,
 Screen-saver function, standby function

Standard accessories



Polarized filters



Remote



Wall mount



Power Adapter

ZN-230 LED VISUAL CHART 23"HD

ECONOMICAL



CE FDA

Both table stand and wall mount are as standard accessories

SPECIFICATIONS

LED monitor.....23-inch TFT Panel, 1920x1080 pixels
 Brightness.....200CD/m²
 Optotype.....E, C(Landolt),Alphabet, Numbers, Children
 Visual unit.....Decimal,Snellen(ft,20/20),Snellen(m,6/6)
 Distance.....1M~6M (in steps of 0.1m)
 Mask type.....Single mask, Horizontal mask, Vertical mask,
 Red/Green mask, White/Black mask,
 Contrast adjustment, Red/Green adjustment
 Special test masks.....Ishihara, Cross cylinder, Astigmatism,
 Cross Grid, Fixation point Test, Red/Green Panel,
 Worth Test, Red/Green Balance, Traffic light,
 Coincidence Test, Stereo Test, Schober Test, ETDRS
 Other functions.....Ramdomization Optotypes Playing function,
 R/G setting function, Video playing function,
 Screen-saver function, standby function



X-22

Optional floor stand for Visual chart,
for screen 19" ~40"

Standard accessories



Wall mount



Remote