

# L3000 Slit Lamp Microscope

## NIPPON SLIT LAMP

### SPECIFICATIONS:

- Microscope: **Stereoscopic**
- Magnification: **10X,16X**
- Eyepiece: **10X**
- Slit Width: **0-12mm**
- Slit Height: **1-12mm**
- Slit Apertures: **φ12、φ9、φ5、φ3、φ1、φ0.2 (mm)**
- Filters: **Heat Absorption, Grey, Redfree, Cobalt Blue**
- Slit Angle: **0°-180°**
- Slit Illumination: **6V20W, Halogen lamp**
- Power source: **A.C.110V or 220V, 50/60Hz**
- Net dimension: **660×285×325mm**
- Net weight: **12Kg**
- Package dimensions: **690×530×235mm**
- Gross weight: **15kg**

