


Cryo-Line



OPTIKON™
Men and Technology

REF 131002

- Simple and automatic system for cryosurgery.
- Wide graphic display.
- User friendly control software.
- Probe temperature control.
- Count down function.
- High precision temperature reading.
- System efficiency test.
- Suitable for retinal, endovitrear, glaucoma and cataract surgery.
- Wide range of cryoprobes.



TECHNICAL SHEET

TECHNICAL DATA

- **GAS TYPE**
N₂O or CO₂ medical grade in non syphon cylinder
- **WORKING PRESSURE**
50/60 bar (725-870psi)
- **MAX PRESSURE**
70 bar (1015psi)

POWER SUPPLY

- **INPUT VOLTAGE**
100/120/220/240V (User Selectable)
- **FREQUENCY**
50/60Hz
- **POWER CONSUMPTION**
80VA

DIMENSIONS

- **HEIGHT**
172mm
- **WIDTH**
393mm

- **DEPTH**
463mm
- **WEIGHT**
10Kg

CLASSIFICATION AND SAFETY STANDARDS

- **COMPLIES**
EN60601-1, EN60601-1-1, EN60601-1-2, 93/42/CEE
- **CLASSIFICATION**
Class I, Type BF (Applied part)
- **PROTECTION AGAINST THE INGRESS OF WATER**
Console IPX1
Footswitch: IPX8

CRYO PROBES



- 133001 Cryo probe for cataract



- 133002 Cryo probe for retina, ball



- 133003 Cryo probe for retina, spatula



- 133004 Cryo probe for retina, "T" shaped



- 133005 Cryo probe for retina, "Bonnet"



- 133006 Cryo probe for glaucoma



- 133007 Cryo probe, endovitrear

MISCELLANEOUS



- 132005 Spare cartridge for gas filter



- 132006 CO₂ cylinder adapter (Chinese standard)



- 132008 N₂O cylinder adapter (ISO5145 standard)



- 132009 CO₂ cylinder adapter (ISO5145 standard)



- 132010 N₂O cylinder adapter (French standard)



- 183001 Trolley for Cryoline with cylinder holder

EUROTECH OPTICAL

For the best Prices &

Quality on Optical Equipment,
Optical Frames & Contact Lenses

PO BOX 13893 Laudium 0037

Pretoria South Africa

Tel: (012) 370 4175

Fax: (012) 370 1557

EUROTECH OPTICAL

For the best Prices &

Quality on Optical Equipment,
Optical Frames & Contact Lenses

PO BOX 13893 Laudium 0037

Pretoria South Africa

Tel: (012) 370 4175

Fax: (012) 370 1557

WORLD DEALERS

EUROPE

Austria
Belgio
Danimarca
Francia
Germania
Grecia
Italia
Olanda
Polonia
Portogallo
Regno Unito
Repubblica Ceca
Russia
Slovacchia
Spagna
Svizzera
Turchia
Ucraina
Ungheria

AFRICA

Algeria
Egitto
Marocco
Sud Africa
Tunisia

NORTH AMERICA

Canada
USA

CENTRAL AMERICA

Cuba
Messico

SOUTH AMERICA

Brasile
Cile
Colombia
Ecuador
Perù

MIDDLE EAST

Arabia Saudita
Bahrain
Georgia
Giordania
Iran
Israele
Kuwait
Libano
Pakistan
Siria
Turchia

ASIA

Bangladesh
Cina
Corea del Sud
Giappone
Hong Kong
India
Indonesia
Tailandia
Taiwan



Distributed by:



TM denotes trademark of Optikon 2000 S.p.A.
All right reserved.

