

LED CHART PROJECTOR TCP-1000

HIGH CONTRAST PROJECTOR

DELIGHT IN SIGHT

Plug and project forever and
ever. No change of lamps.
No maintenance.



- Maintenance-free LED light
- Brighter and clearer chart display
- Sophisticated design
- Quick chart selection
- Programmable display sequence
- Cyrillic chart disc available





 **TOMEY**
TECHNOLOGY AND VISION

THE TOMEY LED CHART PROJECTOR TCP-1000



QUALITY IN DETAIL

The Tomey LED Chart Projector **TCP-1000** is designed for durability and efficient operation. Thanks to the cold LED light source the lamp is maintenance-free due to its up to 500.000 hours lifetime. The intuitive remote control allows you to program three customised display sequences. The **TCP-1000** will be delivered including remote control, projection screen and base for immediate use. You can also easily mount the system on your refraction unit or onto the wall.



SPECIFICATIONS

MEASUREMENT RANGES

Auto-off function	Adjustable 5 / 10 / 15 minutes
Working distance	2.9 to 6.1 m
Projection magnification	30 x
Projection tilt	Approx. $\pm 15^\circ$
Background illumination	200 cd/m ² at 5 m (adjustable)
Binocular vision test	Red-green / polarising
Available charts	36
Masks	Single, horizontal, vertical, Red-green / polarising
Optotypes (standard)	Numbers, letters, pictures, Snellen E, Landolt
Optional disk	Cyrillic letters

DIMENSIONS & ELECTRIC REQUIREMENTS

Dimensions WDH	191 x 310 x 200 mm
Weight	4 kg
Lamp	3V LED (373 to 801 nm)
Input	AC 115 V to 230 V
Frequency	50/60 Hz
Power consumption	11 VA



Remote control

