

Claw Chopper




 1.5mm D30-1446	 1.0mm D30-1447	 0.7mm D30-2327	 0.9mm D30-2329
 0.7mm D30-2332	 1.0mm D30-2333	 0.7mm D30-2370	 1.0mm D30-2375
 1.6mm D30-2377	 1.6mm D30-2386	 0.8mm D30-2392	 0.8mm D30-2393
 0.8mm D30-2395	 0.8mm D30-2396		



D30-2390

Choppers: Claw Chopper



Claw Chopper

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

 [Click here to view the pdf brochure](#)

Claw Chopper

1.5mm claw shape

D30-1446

Green nucleus divider cutting 45 degree to access, 65 degree angle shaft, tip to angle length 10mm

D30-1447

Ball tip 0.7mm inferior wedge

D30-2327

Naghara Phaco Chopper Rounded, blunt tip 0.9mm inferior wedge

D30-2329

DT30-2329 titanium