


HRR Standard Pseudoisochromatic Test



HRR Standard Pseudoisochromatic Test

Rating: Not Rated Yet  
[Ask a question about this product](#)

Description

[Click here to view the pdf catalogue](#)

Quickly screens for normal color vision and congenital and acquired color vision defects, classifying to type—protan/deutan (red/green) and tritan (blue/yellow)—and extent of defect. This test is ideal for pediatric testing, because it uses simple shapes (O, X, ?) that are easily recognized by children. Special non-glare laminate keeps colors from becoming distorted and prevents fingerprint smudges. Reversible occluding glasses are included.

Unlaminated . . . . .	.730006
Laminated . . . . .	.730005