

## Medop SL-700-2H



Halogen 2 Steps Slit Lamp

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

[Click here to view the pdf brochure](#)

### Halogen 2 Steps

- The SL-700 series slit lamp has 3 different models, 2-step, 3-step and 5-step.
- The resolution is improved to 180 lp/mm, it can give you a clean and clear image.
- The illuminance is greatly improved to over 280Klux. It meets the demand of Fundus diagnosis. Especial when these two symptoms happen to the patients: crystalline lens opacities or vitreous opacities.
- Comfortable mechanical operating feeling.
- The range of refractive degree of eyepiece is improved. Now is from +6.00D to -8.00D. When tuning the refractive degree of eyepiece, the cover of eyepiece does not turn along with it, so that avoid the inconvenience for user.
- The distance of left and right optics is 23mm, which is compatible to almost all the brands of beam splitters in the market.
- It's compatible with the Goldmann applanation tonometer BQ type.

- The compatibility to different brands of laser adaptor is improved.