

P.L. Slit Lamp HS-02



Two Step Magnification

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Slit Lamp H-Series

[Click here to view the pdf catalogue](#)

Binocular Microscope

Eye Pieces: 10x, 15x

Objective: 1x, 1.6x

Magnification: 10x, 16x and 24x

Illumination Unit

Slit Image Rotation: 0 to 180°

Tilting Illumination: Up to 20°

Filter Disc.: Cobalt Filter, Green Filter, Yellow Filter,
Natural density and open aperture

Slit Diaphragm Disc.: Six apertures of 12, 9, 7, 3 and 0.2 mm
and a wedge shaped diaphragm to infinitely variable slit lengths
Halogen Lamp: 12 Volts, 50 Watts

Standard Accessories

Replacing Lamp - one piece for illumination unit.
Replacing Mirror
Plastic breath shield
Testing Rod
Dust Cover and Fuse