

TC-700 OPHTHALMIC UNIT & CHAIR



Non-recliner with double top table.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

 [Click here to view the pdf catalogue](#)

FEATURES & SPECIFICATIONS

- Small and exquisite , can be used for small room
- Easy operation with a well positioned control panel
- Low voltage output for ophthalmoscope and Retinoscope voltage selectable, 2.SV, 3V, 3.SV, 4V, 6V, 12V
- "Keller" type and "Welch Allyn" type charger are for optional
- Projector mounting plate
- Table and case using composite material