

Eurotech VS-2000I



Optovision Chair and Stand with double top table.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Features

Big sliding table to accommodate 2 instruments

3 rechargeable wells available

Chair can be reclined to 165°

Specifications

Base size.....39inch x 32inch (1000mm x 800mm)

Tabletop size.....33inch x 16inch (840mm x 405mm)

Sliding table to accommodate two instruments

Auto "power off" for table instruments.

Table instrument voltage selectable: 6V and 12V)

Auxiliary 110V/220V output under table

Three rechargeable wells

Counterbalanced Refractor Arm

(a) Adjustable Load Range: 8lbs to 20 lbs. (3.6kgs to 9kgs)
(b) Vertical Adjustment Range: 12inch (405mm)
Table rotate angle.....90°
Chair rotate..... ±90°
Recline angle of chair.....165°
Chair height.....20inch ? 25inch (500mm ? 640mm)
Adjustable range of headrest.....4inch (100mm)
White LED lamp..... 10 W, brightness is adjustable
Voltage.....110/220V. (Switch voltage freely.)
Gross weight.....485lbs (220kgs)